

APR Tables
CHIRO - 400FS Techniques
48" 20/20 Digital Panel

SID:	72"
Focal Spot:	Large

AP														
cm	16	18	20	22	24	26	28	30	32	34	36	38	40	42
kVp	80	80	80	80	80	85	85	90	90	90	90	90	95	95
mAs	20	30	40	50	60	70	80	95	110	130	150	180	180	200

Lateral														
cm	18	20	22	24	26	28	30	32	34	36	38	40	42	44
kVp	80	80	80	80	82	82	85	85	90	93	95	95	98	98
mAs	50	60	75	80	85	95	100	120	130	140	160	180	180	200

Lateral Cervical / Thoracic														
cm	18	20	22	24	26	28	30	32	34	36	38	40	42	44
kVp	80	80	80	80	80	82	85	85	90	90	93	93	95	95
mAs	40	50	60	75	85	95	110	120	120	140	150	170	180	200

Lateral Thoracic / Lumbar														
cm	18	20	22	24	26	28	30	32	34	36	38	40	42	44
kVp	80	80	80	83	85	85	85	90	90	90	95	95	95	100
mAs	50	60	70	80	95	100	110	120	140	160	160	180	200	200